



Split Gate

Specifications

Assist with the Calving Process Artificial Insemination Pen Sick Pen Can be Set Up Inside a Building for a More Permanent Installation Simple to Operate Low-Cost Investment Built Strong to Last for Many Years Safety for Operator & Animal is Top Priority 7-Bar (69" tall) Side Panel Height Built of 1 1/4 " & 2" Square Tubing Split Squeeze Gate to Allow Nursing or for C-Sections, Sheet Metal on Bottom Half of Gate Stanchion Headgate with a Rope Close & Emergency Quick Release Standard Size: 10' x 10' 7" Bar Pen w/ Stanchion (69" tall)

T.W.I. Calving Pen

Simple, Safe, Portable & Built to Last



Quick Release Adjustable Headgate



Calving Pen Latch

To Contact Your Local Dealer! call: 800.252.0847 online: www.titanwestinc.com